



THE WORLD'S FIRST DNA/RNA DECONTAMINATION TABLET

DNA/RNA Decontamination. Reimagined.

One tablet transforms 250 ml of distilled water
into a fresh, standardized working solution.



1 tablet = 250 ml


WORKING SOLUTION

< half a minute

CONTACT TIME

24 months

TABLET SHELF LIFE



INTACT DNA/RNA

FRAGMENTED

One system. Three formats.

From single laboratories to centralized supply

10 TABLETS



Jar 10 x 250 ml DN-X

2.5 L working solution

CODE 90001730

250 TABLETS



Box 250 x 250 ml DN-X

62.5 L working solution

CODE 90001740

500 TABLETS



Box 500 x 250 ml DN-X

125 L working solution

CODE 90001741

FRESHLY PREPARED · NON-CORROSIVE · BIODEGRADABLE · COMPACT SHIPPING

DN-X™ BY ALLVATER BIOSOLUTIONS

For laboratories, distributors, OEM and private label

Where to order?

contact:

info@av-biosolutions.de

www.dn-x.de

www.av-biosolutions.de